EXPERIENCE THE PERFORMANCE

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


Group difference (95% confidence interval):
-2.78 (-3.44 to -2.13) percent-seconds;
data reflect least squares mean ± standard error.

Group difference (95% confidence interval):
-6.12 (-9.82 to -2.43) mL; data reflect least squares mean ± standard error.

Group difference (95% confidence interval):
-15.78 (-26.49 to -5.07) seconds; data reflect least squares mean ± standard error.

Centurion® Configuration:
Centurion® Vision System,
45° Balanced Tip with Intrepid® Ultra Sleeve

INFINITI® Configuration:
INFINITI® Vision System, 45° Mini Flared
Kelman with Ultra Sleeve

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


The Centurion® Vision System handpiece featuring the INTREPID® Balanced Tip and INTREPID® Sleeve

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BALANCED Tip
technology

¹As compared to the INFINITI® Vision System, bottle gravity system.

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*Trademarks are the property of their respective owners.

1. Zacharias J. Thermal characterization of phacoemulsification probes operated in axial and torsional modes. *J Cataract Refract Surg.* 2015;41(1):208-216.

2. Han YK, Miller KM. Heat production: longitudinal versus torsional phacoemulsification. *J Cataract Refract Surg.* 2009;35(10):1799-1805.

3. Christakis PG, Braga-Mele RM. Intraoperative performance and postoperative outcome comparison of longitudinal, torsional and transversal phacoemulsification machines. *J Cataract Refract Surg.* 2012;38(2):234-241.

4. Sharif-Kashani P, Fanneq D, Injev V. Comparison of occlusion break responses and vacuum rise times of phacoemulsification systems. *BMC Ophthalmol.* 2014;14:96.

5. Nicoli M, Miller K, Dimalanta R, Loke D; Jules Stein Eye Institute, UCLA. IOP Stability Measurement and Comparison Between Gravity-Fed and Actively Controlled Phacoemulsification Systems. 2014.

6. Solomon K, Lorente R, Cionni R, Fanneq D. Prospective, randomized clinical study using a new phaco system with intraocular system target pressure control. ASCRS-ASOA Symposium and Congress; April 25-29, 2014; Boston, USA.

7. Alcon data on file.

8. Rowen S. Microincision cataract surgery. Presented January 2014; Park City, UT.

9. Boukhny M, Sorensen G, Gordon R. A novel phacoemulsification system utilizing feedback based IOP target control. ASCRS Presentation. 2014.

10. Zacharias J. Comparative thermal characterization of phacoemulsification probes operated in elliptical, torsional and longitudinal ultrasound modalities. ASCRS-ASOA Symposium and Congress; April 25-29, 2014; Boston, USA.

11. Zacharias J. Comparative motion profile characterization of the mini flared and balanced phacoemulsification tips. ESCRS; September 5-9, 2015; Barcelona, Spain.

12. Vasavada AR, et al. Comparison of torsional and microburst longitudinal phacoemulsification: a prospective, randomized, masked clinical trial. *Ophthalmic Surg Lasers Imaging.* 2010;41(1):109-114.

13. Allen D, Habib M, Steel D. Final incision size after implantation of a hydrophobic acrylic aspheric intraocular lens: new motorized injector versus standard manual injector. *J Cataract Refract Surg.* 2012;38(2):249-255.

14. Johansson C. Comparison of Motorized IOL Insertion to Traditional Manual IOL Delivery. ASCRS; March 25-29, 2011; San Diego, USA.

Alcon[®]

a Novartis company